


# SUREHOOK


## Vitrification Straws



# SUREHOOK

## Vitrification Straws


The SUREHOOK Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.




Hook tip



Flat surface wide  
body for writing



Transparent  
Sleeve



Ring marking for easy  
visualization during insertion

| Colour    | White | Red   | Green | Blue  | Yellow | Qty. / Box |
|-----------|-------|-------|-------|-------|--------|------------|
| Order No. | SUH-W | SUH-R | SUH-G | SUH-B | SUH-Y  | 10         |

\*All products may not be available in all countries. Please contact our sales department for details.



### US Sales

✉ sales.usa@allwinmedical.com  
☎ +1 714-572-1709

### India Sales

✉ sales.india@allwinmedical.com  
☎ +91 99870-84563

### Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com  
☎ +91 95949-81581