



## Injection Needle Metallic



The Injection Needle Metallic is used for cystoscopic injection into urethra, bladder neck & the bladder wall.

- Bevel tip designed for atraumatic placement and ease of advancement
- Needle luer lock
- Needle constructed of stainless steel
- Circular marking provided for accurate positioning and placement of needle
- Straight marking to show needle bevel position
- Offered in 23G



Order No.	Gauge	Length (cm)	Qty. / Box
INJ-M23/35	23	35	5

\*All products may not be available in all countries. Please contact our sales department for details.