



Krait Ureteral Catheter - Closed Tip



The Krait Ureteral Catheter - Closed Tip is used for drainage of urine from the ureter and the delivery of contrast media.

- Closed tip allows placement over a stylet for negotiating a tortuous ureter
- Radiopaque catheter to enable accurate positioning
- Centimetre marking provided for precise advancement and placement of the catheter
- Supplied with stainless steel stylet & ureteral catheter adaptor
- Soft polyurethane material provides superior biocompatibility and drainage. It also ensures maximum patient comfort



Order No.	Fr	Length (cm)	Qty. / Box
K-UC-R3	3	70	10
K-UC-R4	4	70	10
K-UC-R5	5	70	10
K-UC-R6	6	70	10
K-UC-R7	7	70	10

*All products may not be available in all countries. Please contact our sales department for details.