



ARMADA Amplatz Renal Dilator Set



The ARMADA Amplatz Renal Dilator Set is used for progressive dilatation of the nephrostomy tract prior to percutaneous kidney stone removal.

- Radiopaque dilators & sheath to enable accurate positioning
- Gradual taper tip designed for atraumatic placement and maneuverability
- Smooth surface assists in easier introduction of catheter
- Kink resistant
- Introducer catheter acts as a guide for dilator 14 Fr to 30 Fr



Order No. (Amplatz Set)	Description	Qty. / Box
ARDS8/30	Amplatz Renal Dilatation Set Includes : 1) 8, 10, 12, 14, 16 Fr - PE Fascial Dilators 2) 18, 20, 22, 24, 26, 28, 30 Fr - PE Renal Dilators 3) 28, 30 Fr Amplatz Sheath - PTFE 4) Introducer Catheter for Dilators 8 Fr, 70 cm	1
Order No. (Amplatz Full Set)	Description	Qty. / Box
ARDS8/30-FS	Amplatz Renal Dilatation Full Set Includes : 1) 8, 10, 12, 14, 16 Fr - PE Fascial Dilators 2) 18, 20, 22, 24, 26, 28, 30 Fr - PE Renal Dilators 3) 24, 26, 28, 30 Fr Amplatz Sheath - PTFE 4) Introducer Catheter for Dilators 8 Fr, 70 cm 5) Introducer Catheter Sheath PTFE 10 Fr, 30 cm	1

*All products may not be available in all countries. Please contact our sales department for details.