


ARMADA Nottingham One Step Dilator



The ARMADA Nottingham One Step Dilator is Used for one step dilatation of the intramural portion and ureter.

- One step dilation for the intramural portion and ureter. Intended for one time use only
- May also be used to inject using female luer lock adaptor at proximal end
- Designed to be used over 0.038" guidewire
- 4cm taper length, overall length is 60cm

Integrally luer lock
molded hub

Radiopaque tubing

Long taper tip



Order No.	Tapered Size Fr	Guidewire	Qty. / Box
NOT6/10	6	0.038"	10
NOT6/12	6	0.038"	10

*All products may not be available in all countries. Please contact our sales department for details.