



ARMADA Renal Dilator with Amplatz Sheath



The ARMADA Renal Dilator with Amplatz Sheath used for nephrostomy tract dilation and the placement of sheath.

- Radiopaque dilators & sheath to enable accurate positioning
- Gradual tapered tip designed for atraumatic placement and maneuverability
- Smooth surface assists in easier introduction of catheter
- Kink resistant
- Introducer catheter acts as a guide for dilator of 14 Fr to 30 Fr



Order No.	Fr	Length (cm)	Qty. / Box
RD-AM12/30	12	30	5
RD-AM14/30	14	30	5
RD-AM16/30	16	30	5
RD-AM18/30	18	30	5
RD-AM20/30	20	30	5
RD-AM22/30	22	30	5
RD-AM24/30	24	30	5
RD-AM26/30	26	30	5
RD-AM28/30	28	30	5
RD-AM30/30	30	30	5

*All products may not be available in all countries. Please contact our sales department for details.