



ARMADA Renal Dilator



The ARMADA Renal Dilator is used for dilatation of renal tract.

- Radiopaque dilator to enable accurate positioning
- Gradual taper tip designed for atraumatic placement and maneuverability
- Smooth surface assists in easier introduction of guidewire
- Wide range of sizes available



Order No.	Fr	Length (cm)	Qty. / Box
RD12	12	30	5
RD14	14	30	5
RD16	16	30	5
RD18	18	30	5
RD20	20	30	5
RD22	22	30	5
RD24	24	30	5
RD26	26	30	5
RD28	28	30	5
RD30	30	30	5

*All products may not be available in all countries. Please contact our sales department for details.