



ARMADA Ureteral Dilators



The ARMADA Ureteral Dilators is used for step by step dilatation of the intramural tunnel and ureter and/or for stone manipulation.

- Designed to be used over a 0.038" Guidewire
- 60 cm length with luer adaptor to permit injection
- Gradual tapered tip facilitates easy insertion and maneuverability of dilator into the ureter



Order No. (Ureteral Dilator Set)	Length (cm)	Description	Qty. / Box
UDS6/16-60	60	Ureteral Dilator Set of (6F, 8F, 10F, 12F, 14F, 16F)	1

Order No. (Individual Ureteral Dilator)	Fr	Length (cm)	Guidewire	Qty. / Box
UD6	6	60	0.038	10
UD8	8	60	0.038	10
UD10	10	60	0.038	10
UD12	12	60	0.038	10
UD14	14	60	0.038	10
UD16	16	60	0.038	10

*All products may not be available in all countries. Please contact our sales department for details.