


Advanced Inflation Device



The Advanced Inflation Device is used with balloon dilation catheter to monitor inflation pressure.

- Ergonomic design facilitates excellent eye to hand co-ordination
- Efficient intra-cam locking mechanism allows to maintain pressure and guarantees fast and automatic release
- Superior design with grained surface and slip free lever ensures maximum grip even with wet hands
- Pressure gauge ensures fatigue free readability
- Upto 30/26 atm analogue pressure range
- Registers immediately without pressure decay or drift
- Easy locking and unlocking



Order No.	Description	Qty. / Box
IDA20	Inflation Device Advanced 20cc	1

*All products may not be available in all countries. Please contact our sales department for details.