



BSST - Biliary Stents Standard PTFE



The BSST - Biliary Stents Standard PTFE are used for drain obstructed biliary ducts.

- Radiopaque PTFE tubing material
- Smooth rounding posterior
- Taper tip for easy insertion
- Various sizes available



Order No.	Fr	Length (cm)	Qty. / Box
BSST7/5	7	5	10
BSST7/7	7	7	10
BSST7/10	7	10	10
BSST7/12	7	12	10
BSST7/15	7	15	10
BSST10/5	10	5	10
BSST10/7	10	7	10
BSST10/10	10	10	10
BSST10/12	10	12	10
BSST10/15	10	15	10

Order No. (Kit)	Fr	Length (cm)	Qty. / Box
BSST7/5-K	7	5	1
BSST7/7-K	7	7	1
BSST7/10-K	7	10	1
BSST7/12-K	7	12	1
BSST7/15-K	7	15	1
BSST10/5-K	10	5	1
BSST10/7-K	10	7	1
BSST10/10-K	10	10	1
BSST10/12-K	10	12	1
BSST10/15-K	10	15	1

*All products may not be available in all countries. Please contact our sales department for details.

Kit includes :

- Biliary Stents
- Majestick - Biliary Stent Introducer
- 0.035" PTFE Guidewire 400cm