



Chiba Needle - 2 Part Bevel Tip



The Chiba Needle - 2 Part Bevel Tip is used for percutaneous cholangiography cytological biopsy and routine percutaneous application.

- Bevel tip
- Stainless steel
- Wide range available



Order No.	Gauge	Length (cm)	Qty. / Box
CH-18-15-2	18 G	15	5
CH-18-20-2	18 G	20	5
CH-20-15-2	20 G	15	5
CH-20-20-2	20 G	20	5
CH-22-15-2	22 G	15	5
CH-22-20-2	22 G	20	5
CH-22-28-2	22 G	28	5

*All products may not be available in all countries. Please contact our sales department for details.