


## Cobra Nitinol Guidewire



The Cobra Nitinol Guidewire is used for percutaneous entry into vessel using the seldinger technique or endoscopic access.

- Kink resistant nitinol shaft
- Black and white stripes for easier visualization
- The guidewire provides an ideal blend of stiffness and flexibility for effective control and navigation through a tortuous anatomy
- Straight flexible tip



Order No.	Diameter	Length (cm)	Qty. / Box
CGS150-25	0.025"	150	5
CGS150-32	0.032"	150	5
CGS150-35	0.035"	150	5
CGS450-25	0.025"	450	1
CGS450-32	0.032"	450	1
CGS450-35	0.035"	450	1

\*All products may not be available in all countries. Please contact our sales department for details.