



ENJFT - Endoscopic Naso Jejunal Feeding Tube



The ENJFT - Endoscopic Naso Jejunal Feeding Tube is used through channel of endoscope to provide short term enteral access for delivery of medication or nutrition to the small bowel through an endoscope.

- Soft, biocompatible polyurethane tubing material
- Compatible with 0.035" guidewire
- 5 Centimetre marking
- Radiopaque tubing for easy visualization
- Taper tip for easy insertion



Order No.	Description	Fr	Guidewire	Length (cm)	Qty. / Box
ENJFT8	Endoscopic Naso Jejunal Feeding Tube	8	0.035"	240	5
ENJFT10	Endoscopic Naso Jejunal Feeding Tube	10	0.035"	240	5

Order No.	Description (Kit)	Fr	Guidewire	Length (cm)	Qty. / Box
ENJFT8-K	Endoscopic Naso Jejunal Feeding Tube Kit	8	0.035"	240	1
ENJFT10-K	Endoscopic Naso Jejunal Feeding Tube Kit	10	0.035"	240	1

*All products may not be available in all countries. Please contact our sales department for details.

Each Kit includes :

- 1 ea Feeding Tube
- 1 ea Guidewire 0.035"
- 1 ea Nasal Tube