


## Ellick Bladder Evacuator



The Ellick Bladder Evacuator is used for bladder wash during transurethral resection of prostate (TURP).

- Reusable
- Adapters connects to storz, wolf and acmi
- Silicone bulb eliminates latex allergy



Order No.	Description	Qty. / Box
EE-G	Ellick Evacuator - Glass	1
EE-GR	Ellick Glass Reservoir	1
EE-CT	Ellick Connector with Tubing	1
EE-SB	Ellick Silicone Bulb	1

\*All products may not be available in all countries. Please contact our sales department for details.