


Fascial Dilators



The Fascial Dilators are used for dilatation of the percutaneous tract over the guidewire.

- Radiopaque tubing made up of polyethylene material to maximize patient comfort
- Female luer lock allows heparin, contrast media and other solutions to be injected using standard luer lock syringe
- Uniform taper minimizes trauma to patient
- Smooth surface assists in easier introduction of guidewire
- Integrally molded hub



Order No. (Fascial Dilator Set)	Length (cm)	Description	Qty. / Box
FDS 6/10	20	Fascial Dilator Set (6 Fr, 7 Fr, 8 Fr, 9 Fr, 10 Fr)	5
FDS 8/16	20	Fascial Dilator Set (8 Fr, 10 Fr, 12 Fr, 14 Fr, 16 Fr)	5
FDS 6/16	20	Fascial Dilator Set (6 Fr, 8 Fr, 10 Fr, 12 Fr, 14 Fr, 16 Fr)	5

Order No. (Individual Fascial Dilator)	Fr	Length (cm)	Guidewire	Qty. / Box
FD6	6	20	0.038	5
FD7	7	20	0.038	5
FD8	8	20	0.038	5
FD9	9	20	0.038	5
FD10	10	20	0.038	5
FD12	12	20	0.038	5
FD14	14	20	0.038	5
FD16	16	20	0.038	5

*All products may not be available in all countries. Please contact our sales department for details.