


FlexiGuide Standard - Straight Tip Guidewire



The FlexiGuide Standard - Straight Tip Guidewire is used for percutaneous entry into vessel using the seldinger technique or for endoscopic access.

- Straight tip fixed core
- PTFE coating allows smooth catheter insertion and stent tracking

PTFE coating

Straight tip

| Order No. | Tip | Diameter | Length (cm) | Qty. / Box |
|-----------|----------------------------------|----------|-------------|------------|
| SFT150-18 | Straight Fixed Core, PTFE Coated | 0.018" | 150 | 5 |
| SFT150-25 | Straight Fixed Core, PTFE Coated | 0.025" | 150 | 5 |
| SFT150-28 | Straight Fixed Core, PTFE Coated | 0.028" | 150 | 5 |
| SFT150-32 | Straight Fixed Core, PTFE Coated | 0.032" | 150 | 5 |
| SFT150-35 | Straight Fixed Core, PTFE Coated | 0.035" | 150 | 5 |
| SFT150-38 | Straight Fixed Core, PTFE Coated | 0.038" | 150 | 5 |

*All products may not be available in all countries. Please contact our sales department for details.