



## GRAB-IT Stainless Steel Stone Basket - Flatwire & Helical



The GRAB-IT Stainless Steel Stone Basket has been designed for endoscopic manipulation, dislodging, engaging and removing renal, ureteral or biliary stones.

- Constructed of stainless steel - 4 wires
- Offered in sizes 3 Fr and 4 Fr to meet all requirements
- Ergonomic handle provides control and convenient manipulation
- Dismantlable handle



Order No. (Flatwire)	Order No. (Helical)	Fr	Length (cm)	Qty. / Box
SF4W/3F-90-F	SB4W/3F-90-H	3	90	1
SF4W/3F-120-F	SB4W/3F-120-H	3	120	1
SF4W/4F-90-F	SB4W/4F-90-H	4	90	1
SF4W/4F-120-F	SB4W/4F-120-H	4	120	1

\*All products may not be available in all countries. Please contact our sales department for details.