


GRASP-IT Three - Prong Grasping Forcep



The GRASP-IT Three - Prong Grasping Forcep is designed for endoscopic manipulation, dislodging, engaging and removal of renal and ureteral stones.

- Made up a stainless steel
- Sharp prong 3 wire
- Flexible shaft for easier manipulation



Order No.	Fr	Length (cm)	Qty. / Box
GF3W/3F-90	3	90	1
GF3W/3F-120	3	120	1
GF3W/4F-90	4	90	1
GF3W/4F-120	4	120	1

*All products may not be available in all countries. Please contact our sales department for details.