


HyGuide Hydrophilic Guidewire



The HyGuide Hydrophilic Guidewire is used for access ureter during routine and difficult cases.

- Exacting tolerance of each wire facilitates easy and reliable instrument delivery
- Hydrophilic coating permits smooth navigation through tortuous anatomy
- Non-kink shaft



Order No.	Tip	Diameter	Coated	Length (cm)	Qty. / Box
HGSF150-35	Straight Fixed Core	0.035"	Hydrophilic	150	5

*All products may not be available in all countries. Please contact our sales department for details.