


Initial Puncture Needle - 2 Part Trocar Tip



The Initial Puncture Needle - 2 Part Trocar Tip is used for percutaneous nephrostomy and all other percutaneous applications where a 0.035" guidewire is to be used.

- Trocar tip
- Sharp stainless steel tip for increased control and easier penetration through the fascia



Order No.	Tip	Gauge	Length (cm)	Qty. / Box
IPN-18-20-2T	2 Part Trocar	18	20	5
IPN-18-30-2T	2 Part Trocar	18	30	5
IPN-18-15-2T	2 Part Trocar	18	15	5

*All products may not be available in all countries. Please contact our sales department for details.