


Krait Ureteral Catheter - Open Tip



Krait Ureteral Catheter - Open Tip is used for drainage of urine from the ureter and the delivery of contrast media.

- Open tip allows placement over a guidewire for negotiating a tortuous ureter
- Radiopaque catheter to enable accurate positioning
- Centimetre marking provided for precise advancement and placement of the catheter
- Supplied with ureteral catheter adaptor
- Soft polyurethane material provides superior biocompatibility and drainage. It also ensures maximum patient comfort



Order No. Open Tip	Fr	Length (cm)	Qty. / Box
K-UC-O3	3	70	10
K-UC-O4	4	70	10
K-UC-O5	5	70	10
K-UC-O6	6	70	10
K-UC-O7	7	70	10

*All products may not be available in all countries. Please contact our sales department for details.