



Maestro Nitinol Tipless Stone Basket



The Maestro Nitinol Tipless Stone Basket has been designed for endoscopically manipulating, dislodging, engaging and removing renal, ureteral or biliary stones.

- Tipless basket design to help with atraumatic retrieval of stones in the ureter and calyces
- Constructed of nitinol, the stone basket offers a kink-resistant and flexible design to help with stone extraction
- Dismantlable handle
- Laser cut nitinol basket designed to maintain basket stability during deployment
- Distal end of our 1.8 Fr stone basket shaft is made thinner for more flexibility



Order No.	Fr	Length (cm)	Qty. / Box
M-TNB 4W/1.8F-120	1.8	120	1

*All products may not be available in all countries. Please contact our sales department for details.