



## Majestick - Biliary Stents Introducer



The Majestick - Biliary Stents Introducer is used for endoscopic biliary placement.

- PTFE material with radiopaque band marker
- 6 Fr guiding catheter
- Polypropylene 10 Fr pusher catheter
- Glides smoothly over guidewire
- Trocar hub with luer lock



Order No.	Application Stent Size	Guiding Catheter	Guiding Catheter Length	Pusher Catheter Size	Pusher Catheter Length	Qty. / Box
BSI10	10 Fr	6 Fr	208 cm	10 Fr	170 cm	5

\*All products may not be available in all countries. Please contact our sales department for details.