


NephroTrack Nephrostomy Track Balloon Dilator



The NephroTrack Nephrostomy Track Balloon Dilator is recommended for dilating the nephrostomy tract.

- Accepts 0.038" guidewire
- 17 atm rated burst pressure
- 10 mm (30 Fr) x 12 cm balloon
- Supplied with or without amplatz sheath
- Option to set with 20cc Inflation device



Order No.	Order No. (with Inflation Device) OD	Balloon Inflated OD	Balloon Length	Shaft Fr	Usable Length	Qty. / Box
NT10/12	NT10/12-IDA	10 mm (30Fr)	12	7	55	1
NTAM10/12	NTAM10/12-IDA	10 mm (30Fr)	12	7	55	1

*All products may not be available in all countries. Please contact our sales department for details.