




PAN Stents



The PAN Stents are used for drainage of obstructed pancreatic ducts.

- Radiopaque PEBAX tubing material
- Rounded end
- Pusher components available separately - P5/170 or P7/170
- Taper tip for easy insertion
- Various types & sizes available



Order No. (Pigtail)	Fr	Length (cm)	Qty. / Box
PSP5/3	5	3	10
PSP5/5	5	5	10
PSP5/7	5	7	10
PSP5/9	5	9	10
PSP5/12	5	12	10
PSP7/3	7	3	10
PSP7/5	7	5	10
PSP7/7	7	7	10
PSP7/9	7	9	10
PSP7/12	7	12	10

Order No. (Curved)	Fr	Length (cm)	Qty. / Box
PSC5/3	5	3	10
PSC5/5	5	5	10
PSC5/7	5	7	10
PSC5/9	5	9	10
PSC5/12	5	12	10
PSC7/3	7	3	10
PSC7/5	7	5	10
PSC7/7	7	7	10
PSC7/9	7	9	10
PSC7/12	7	12	10

Order No. (Pigtail) (Kit)	Fr	Length (cm)	Qty. / Box
PSP5/3-K	5	3	1
PSP5/5-K	5	5	1
PSP5/7-K	5	7	1
PSP5/9-K	5	9	1
PSP5/12-K	5	12	1
PSP7/3-K	7	3	1
PSP7/5-K	7	5	1
PSP7/7-K	7	7	1
PSP7/9-K	7	9	1
PSP7/12-K	7	12	1

Order No. (Curved) (Kit)	Fr	Length (cm)	Qty. / Box
PSC5/3-K	5	3	1
PSC5/5-K	5	5	1
PSC5/7-K	5	7	1
PSC5/9-K	5	9	1
PSC5/12-K	5	12	1
PSC7/3-K	7	3	1
PSC7/5-K	7	5	1
PSC7/7-K	7	7	1
PSC7/9-K	7	9	1
PSC7/12-K	7	12	1

*All products may not be available in all countries. Please contact our sales department for details.

Kit includes :

- Pancreatic Stents
- 0.035" PTFE Guidewire 400cm
- Pusher