



R+ Pediatric Ureteral Stents (Open Tip, Closed Tip)



The R+ Pediatric Ureteral Stent is intended to temporarily facilitate drainage of urine from kidney to bladder by trained physician.

- Low encrustation tendency to facilitate easy removal of the stent
- Soft polyurethane material ensures high patient comfort
- Sizes to fit pediatric needs



Order No. (Open Tip)	Order No. (Closed Tip)	Fr	Length (cm)	Guidewire	Qty. / Box
POS3/10-R+	PS3/10-R+	3	10	0.018"	10
POS3/12-R+	PS3/12-R+	3	12	0.018"	10
POS3/14-R+	PS3/14-R+	3	14	0.018"	10
POS3/16-R+	PS3/16-R+	3	16	0.018"	10
POS4/10-R+	PS4/10-R+	4	10	0.025"	10
POS4/12-R+	PS4/12-R+	4	12	0.025"	10
POS4/14-R+	PS4/14-R+	4	14	0.025"	10
POS4/16-R+	PS4/16-R+	4	16	0.025"	10

*All products may not be available in all countries. Please contact our sales department for details. Contact us for other sizes

* Combinations :

- With Guidewire : Add Suffix - G
- With Suture : Add Suffix - S
- With Clamp : Add Suffix - C
- With Hydrophilic Coating : Add Suffix - H