



R+ Renal Transplant Stent (Open Tip, Closed Tip)



The R+ Renal Transplant Stent is intended to temporarily facilitate drainage of urine from kidney to bladder by trained physician.

- Soft polyurethane material ensures patient comfort
- Low encrustation tendency to facilitate easy removal of the stent
- Made of biocompatible polyurethane
- Sizes to fit renal transplant needs
- Marking to aid in stent advancement and precise placement



Order No. (Open Taper Tip)	Order No. (Closed Taper Tip)	Fr	Length (cm)	Guidewire	Qty. / Box
RTS-O6/12-R+	RTS-C6/12-R+	6	12	0.038"	10
RTS-O6/14-R+	RTS-C6/14-R+	6	14	0.038"	10
RTS-O6/16-R+	RTS-C6/16-R+	6	16	0.038"	10
RTS-O7/12-R+	RTS-C7/12-R+	7	12	0.038"	10
RTS-O7/14-R+	RTS-C7/14-R+	7	14	0.038"	10
RTS-O7/16-R+	RTS-C7/16-R+	7	16	0.038"	10

*All products may not be available in all countries. Please contact our sales department for details. Contact us for other sizes

*** Combinations :**

- With Guidewire : Add Suffix - G
- With Suture : Add Suffix - S
- With Clamp : Add Suffix - C
- With Hydrophilic Coating : Add Suffix - H