



R+ Ureteral Stents (Open Tip, Closed Tip and Intraoperative)



The R+ Ureteral Stents are intended to temporarily facilitate the drainage of urine from kidney to bladder by a trained physician.

- Low encrustation tendency to facilitate easy removal of the stent
- Soft polyurethane material ensures high patient comfort
- Wide lumen for superior drainage
- Radiopaque stent to enable accurate positioning
- Centimetre marking to aid in stent advancement and precise placement
- 4.7 Fr stent with a long pusher for use via ureterscope



Order No. (Open Taper Tip)	Order No. (Closed Taper Tip)	Order No. Intraoperative (Both Ends Closed Taper)	Fr	Length (cm)	Guidewire	Qty. / Box
SOT4.7/16-R+	SS4.7/16-R+	SSI4.7/16-R+	4.7	16	0.035"	10
SOT4.7/24-R+	SS4.7/24-R+	SSI4.7/24-R+	4.7	24	0.035"	10
SOT4.7/26-R+	SS4.7/26-R+	SSI4.7/26-R+	4.7	26	0.035"	10
SOT4.7/28-R+	SS4.7/28-R+	SSI4.7/28-R+	4.7	28	0.035"	10
SOT4.7/30-R+	SS4.7/30-R+	SSI4.7/30-R+	4.7	30	0.035"	10
SOT5/16-R+	SS5/16-R+	SSI5/16-R+	5	16	0.035"	10
SOT5/24-R+	SS5/24-R+	SSI5/24-R+	5	24	0.035"	10
SOT5/26-R+	SS5/26-R+	SSI5/26-R+	5	26	0.035"	10
SOT5/28-R+	SS5/28-R+	SSI5/28-R+	5	28	0.035"	10
SOT5/30-R+	SS5/30-R+	SSI5/30-R+	5	30	0.035"	10
SOT6/16-R+	SS6/16-R+	SSI6/16-R+	6	16	0.038"	10
SOT6/24-R+	SS6/24-R+	SSI6/24-R+	6	24	0.038"	10
SOT6/26-R+	SS6/26-R+	SSI6/26-R+	6	26	0.038"	10
SOT6/28-R+	SS6/28-R+	SSI6/28-R+	6	28	0.038"	10
SOT6/30-R+	SS6/30-R+	SSI6/30-R+	6	30	0.038"	10
SOT7/16-R+	SS7/16-R+	SSI7/16-R+	7	16	0.038"	10
SOT7/24-R+	SS7/24-R+	SSI7/24-R+	7	24	0.038"	10
SOT7/26-R+	SS7/26-R+	SSI7/26-R+	7	26	0.038"	10
SOT7/28-R+	SS7/28-R+	SSI7/28-R+	7	28	0.038"	10
SOT7/30-R+	SS7/30-R+	SSI7/30-R+	7	30	0.038"	10

*All products may not be available in all countries. Please contact our sales department for details. Contact us for other sizes

* Combinations :

- With Guidewire : Add Suffix - G
- With Suture : Add Suffix - S
- With Clamp : Add Suffix - C
- With Hydrophilic Coating : Add Suffix - H