


R+ Urinary Diversion Stents (Mono-J) (Open Tip, Closed Tip)



The R+ Urinary Diversion Stents (Mono-J) are used for intraoperative placement in the ureter during percutaneous, endoscopic or operative procedures.

- Available in open and closed taper configurations
- With enhanced radiopacity for patient comfort



Order No. (Open Taper Tip)	Order No. (Closed Taper Tip)	Fr	Length (cm)	GW	Qty. / Box
MJO4.7/70-R+	MJC4.7/70-R+	4.7	70	0.035"	10
MJO6/70-R+	MJC6/70-R+	6	70	0.038"	10
MJO7/70-R+	MJC7/70-R+	7	70	0.038"	10

*All products may not be available in all countries. Please contact our sales department for details.

* **Combinations** : With Guidewire : Add Suffix - G