



SCN - Sclerotherapy Needle



The SCN - Sclerotherapy Needle is used for injection into and around bleeding vessels including varices.

- PTFE tubing for frictionless maneuvering
- Smooth bevel tip to facilitate injection
- 2-position handle ensures safety
- Handle clicks on needle extension
- Wide range of sizes available



Order No.	Fr	Length (cm)	Gauge	Qty. / Box
SCN7/160/21	7	160	21	5
SCN7/160/23	7	160	23	5
SCN7/160/25	7	160	25	5
SCN7/240/21	7	240	21	5
SCN7/240/23	7	240	23	5
SCN7/240/25	7	240	25	5

*All products may not be available in all countries. Please contact our sales department for details.