



Super Strong Tumor Stents



The Super Strong Tumor Stent is used for temporary as well as for palliative treatments of internal urinary drainage when tumors are present in ureteral tract.

- Soft kidney and bladder loops provided for patient comfort
- Rigid shaft for maximum compression resistance
- Drainage holes only on the loops
- Radiopaque platinum/iridium (Pt/Ir) marker bands are provided to guide and place catheters during medical procedures
- Compatible with 0.035" guidewire



Order No.	Fr	Length (cm)	Guidewire	Qty. / Box
TS6/24-PA	6	24	0.035"	5
TS6/26-PA	6	26	0.035"	5
TS6/28-PA	6	28	0.035"	5
TS6/30-PA	6	30	0.035"	5

*All products may not be available in all countries. Please contact our sales department for details.

* **Combinations** : With Guidewire : Add Suffix - G • Each stent is supplied with pusher