



Supremo Nitinol Hydrophilic Guidewire



The Supremo Nitinol Hydrophilic Guidewire is used to access ureter during routine and difficult cases.

- Nitinol non-kink core
- Hydrophilic coated
- In straight and angled design
- Super maneuverability



Order No.	Tip	Diameter	Coated	Length (cm)	Qty. / Box
SHSG 150-25	Straight	0.025"	Hydrophilic	150	5
SHSG 150-32	Straight	0.032"	Hydrophilic	150	5
SHSG 150-35	Straight	0.035"	Hydrophilic	150	5
SHAG 150-25	Angled	0.025"	Hydrophilic	150	5
SHAG 150-32	Angled	0.032"	Hydrophilic	150	5
SHAG 150-35	Angled	0.035"	Hydrophilic	150	5

*All products may not be available in all countries. Please contact our sales department for details.