



## U-Flex DUO Ureteral Access Sheath Dual Lumen



The U-Flex DUO Ureteral Access Sheath Dual Lumen is used to establish two continuous working channels thereby reducing ureteral trauma during multiple instrument exchange and endoscopic urological procedures.

- Dual lumen catheter with 2 working channels eliminates multiple catheterization and provides a second working channel for injection
- Kink resistant stainless steel coil reinforced sheath with excellent stiffness aids in navigation via a tortuous anatomy
- Hydrophilic coating for easy insertion of the catheter
- Spontaneous locking mechanism holds the sheath and catheter together in a fixed position



Order No.	OD	ID	Length (cm)	Qty. / Box
UFDL14/12-36	14 Fr	12 Fr	36	1
UFDL14/12-46	14 Fr	12 Fr	46	1

\*All products may not be available in all countries. Please contact our sales department for details.