



U-Flex Ureteral Access Sheath



The U-Flex Ureteral Access Sheath is used to establish a continuous working channel thereby reducing ureteral trauma during multiple instrument exchange.

- Kink resistant sheath with braided design aids in navigation through a tortuous anatomy
- Radiopaque tubing to aids in the positioning of the tip of the access sheath under fluroscopic guidance
- Long tapered dilator for gradual dilation and less trauma to patient
- Hydrophilic coating for easier insertion of catheter
- Dual latch mechanism holds the sheath & catheter together in a fixed position
- Offered in different sizes on request



Order No.	OD	ID	Length (cm)	Qty. / Box
UF14/12-24	14 Fr	12 Fr	24	1
UF14/12-38	14 Fr	12 Fr	38	1
UF14/12-45	14 Fr	12 Fr	45	1
UF14/12-54	14 Fr	12 Fr	54	1
UF12.7/10.7-38	12.7 Fr	10.7 Fr	38	1
UF12.7/10.7-45	12.7 Fr	10.7 Fr	45	1
UF12/10-35	12 Fr	10 Fr	35	1
UF12/10-45	12 Fr	10 Fr	45	1
UF12/10-55	12 Fr	10 Fr	55	1
UF11.5/9.5-35	11.5 Fr	9.5 Fr	35	1
UF11.5/9.5-45	11.5 Fr	9.5 Fr	45	1
UF11.5/9.5-55	11.5 Fr	9.5 Fr	55	1

*All products may not be available in all countries. Please contact our sales department for details.