



## U-Perc Nitinol Stone Basket



The U-Perc Nitinol Stone Basket is used for stone extraction during percutaneous nephrolithotomy (PCNL.)

- Tipless basket design to help with atraumatic retrieval of stones in the ureter and calyces
- Light weight handle for easy maneuverability
- Stainless steel shaft
- Constructed of nitinol, the stone basket offers a kink-resistant and flexible design to help with stone extraction



Order No.	Fr	Length (cm)	Qty. / Box
UP8/38	8	38	1

\*All products may not be available in all countries. Please contact our sales department for details.