



## Ultima Percutaneous Nephrostomy Drainage Catheter



The Ultima Percutaneous Nephrostomy Drainage Catheter is used for Ultima percutaneous placement of a pigtail catheter in the pelvis for nephrostomy drainage.

- A convenient complete set for percutaneous placement of drainage catheter - ultima pigtail design
- Material provides superior biocompatibility and drainage
- Radiopaque tubing for precise placement
- Locking loop design for better anchoring
- Firm polyurethane



Order No.	Order No. (with UBC)	Fr	Length (cm)	Qty. / Box
U-PCN-P8/30	U-PCN-P8/30-U	8	30	10
U-PCN-P9/30	U-PCN-P9/30-U	9	30	10
U-PCN-P10/30	U-PCN-P10/30-U	10	30	10
U-PCN-P12/30	U-PCN-P12/30-U	12	30	10
U-PCN-P14/30	U-PCN-P14/30-U	14	30	10

\*All products may not be available in all countries. Please contact our sales department for details.

### Set includes :

- Ultima Percutaneous Nephrostomy Drainage Catheter 30cm
- Urine Bag Connector
- Stopcock

### \*Combinations :

- Add Suffix - L
- With Locking Loop System : Add Suffix : LS
- With Hydrophilic Coating : Add Suffix : H