



Ultima Percutaneous Nephrostomy Drainage Set



The Ultima Percutaneous Nephrostomy Drainage Set - Pigtail is used for percutaneous placement of a pigtail catheter in the pelvis for nephrostomy drainage.

- Pigtail design helps to secure the position of the catheter
- Locking loop design available for better anchoring
- Radiopaque catheter to enable accurate positioning
- Firm polyurethane material provides superior biocompatibility and drainage. It also ensures maximum patient comfort



Order No.	Fr	With Fascial Dilators	Qty. / Box
U-PDS-P8	8	7, 8, 9	1
U-PDS-P9	9	7, 8, 9, 10	1
U-PDS-P10	10	6, 8, 10, 12	1
U-PDS-P12	12	6, 8, 10, 12, 14	1
U-PDS-P14	14	6, 8, 10, 12, 14, 16	1

*All products may not be available in all countries. Please contact our sales department for details.

Set includes :

- Ultima Drainage Catheter
- Initial Puncture Needle - 3 Part
- J Tip fixed core PTFE Coated 0.038" Guidewire
- Fascial Dilators
- Urine Bag Connector
- Scalpel
- Catheter Fixing Disk
- Stopcock

*Combinations :

- With Locking Loop System : Add Suffix : LS
- With Lunderquist Guidewire instead of PTFE Guidewire : Add Suffix - L
- With Hydrophilic Coating : Add Suffix : H