


## UreTrack Ureteral Balloon Dilator



The UreTrack Ureteral Balloon Dilator is recommended for the dilatation of the urinary tract.

- Accepts 0.035" guidewire
- 5 Fr shaft
- 18 atm rated burst pressure
- 5 mm (15 Fr) balloon
- Supplied with stopcock
- Option of set with 20cc inflation device



Order No.	Order No. (with inflation device)	Balloon inflation OD	Balloon length	Shaft Fr	Length (cm)	Qty. / Box
UBD5/4	UBD5/4-IDA	5mm (15Fr)	4	5	70	1
UBD5/10	UBD5/10-IDA	5mm (15Fr)	10	5	70	1

\*All products may not be available in all countries. Please contact our sales department for details.