



allDrain Percutaneous Nephrostomy Catheter - Malecot



The allDrain Percutaneous Nephrostomy Catheter - Malecot is used for percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage.

- Malecot design helps to secure the position of the catheter
- Made of biocompatible polyurethane to maximize patient comfort
- Radiopaque catheter enables accurate positioning



Order No.	Order No. (with Urine Bag Connector)	Fr	Length (cm)	Qty. / Box
PCN-M10	PCN-M10-U	10	30	5
PCN-M12	PCN-M12-U	12	30	5
PCN-M14	PCN-M14-U	14	30	5
PCN-M16	PCN-M16-U	16	30	5

*All products may not be available in all countries. Please contact our sales department for details.

Set includes :

- Percutaneous Nephrostomy Catheter - Malecot 30cm
- Urine Bag Connector
- Stopcock