



allDrain Percutaneous Nephrostomy Drainage Set-Malecot



The allDrain Percutaneous Nephrostomy Drainage Set - Malecot is used for percutaneous placement of a malecot catheter in renal pelvis for nephrostomy drainage.

- Soft polyurethane material provides superior biocompatibility and drainage. It also ensures maximum patient comfort
- Malecot design helps to secure the position of the catheter and promotes the drainage of dense fluid with high viscosity
- Radiopaque catheter to enable accurate positioning
- 30 cm length
- With flexible pusher



Order No.	Fr	With Fascial Dilators	Qty. / Box
PDS-M10	10	6, 8, 10, 12	1
PDS-M12	12	6, 8, 10, 12, 14	1
PDS-M14	14	6, 8, 10, 12, 14, 16	1
PDS-M16	16	6, 8, 10, 12, 14, 16, 18	1

*All products may not be available in all countries. Please contact our sales department for details.

Set includes :

- Malecot Drainage Catheter 30 cm
- Initial Puncture Needle - 3 Part
- J Tip fixed core PTFE Coated 0.038" Guidewire
- Fascial Dilators
- Urine Bag Connector
- Scalpel
- Catheter Fixing Disc
- Stopcock

*Combinations :

- With lunderquist Guidewire instead of PTFE Guidewire : Add Suffix - L
- With Hydrophilic Coating : Add Suffix : H