



## allDrain Re-Entry Malecot Catheter



The allDrain Re-Entry Malecot Catheter is used for internal as well as external drainage with combination of nephrostomy catheter with a stent.

- Tip design allows guidewire to pass into the ureter
- Polyurethane material provides superior biocompatibility and drainage. It also ensures maximum patient comfort
- Radiopaque catheter enables accurate positioning
- With pusher, stopcock and urine bag connector
- Provided with extension tube for internal drainage



Order No.	Order No. (with Urine Bag Connector)	Fr	Length (cm)	Qty. / Box
RM10	RM10-U	10	30	1
RM12	RM12-U	12	30	1
RM14	RM14-U	14	30	1
RM16	RM16-U	16	30	1
RM18	RM18-U	18	30	1

\*All products may not be available in all countries. Please contact our sales department for details.

**Set includes :**

- Re-Entry Malecot Catheter
- Urine Bag Connector
- Stopcock