



allSoft Multiloop Ureteral Stents (Open Tip, Closed Tip)



The allSoft Multiloop Ureteral Stent is intended to temporarily facilitate drainage of urine from kidney to bladder by trained physician.

- Open and closed tip for both end multiloop
- Open and closed tip for bladder end multiloop
- Universally serves ureteral lengths from 22 cm to 30 cm



Order No. (Open Taper Tip) (Bladder End Multiloop)	Order No. (Closed Taper Tip) (Bladder End Multiloop)	Order No. (Open Tip) (Both End Multiloop)	Order No. (Closed Tip) (Both End Multiloop)	Fr	Length (cm)	Guidewire	Qty. / Box
SOT4.7/ML	SS4.7/ML	SOT4.7/ML-ML	SS4.7/ML-ML	4.7	22-30	0.035"	10
SOT5/ML	SS5/ML	SOT5/ML-ML	SS5/ ML-ML	5	22-30	0.035"	10
SOT6/ML	SS6/ML	SOT6/ML-ML	SS6/ML-ML	6	22-30	0.038"	10
SOT7/ML	SS7/ML	SOT7/ML-ML	SS7/ML-ML	7	22-30	0.038"	10

*All products may not be available in all countries. Please contact our sales department for details. Contact us for other sizes

*** Combinations :**

- With Guidewire : Add Suffix - G
- With Suture : Add Suffix - S
- With Clamp : Add Suffix - C
- With Hydrophilic Coating : Add Suffix - H