



allSoft Tumor Stents



The Tumor Stents are used for temporary as well as for palliative treatments of internal urinary drainage when tumors are present in ureteral tract.

- Full shaft with extended diameter
- Made of biocompatible polyurethane to maximize patient comfort
- Soft and flexible material is used for patient comfort allowing prolonged indwelling period
- Low encrustation tendency to facilitate easy removal of the stent
- Radiopaque stent to enable accurate positioning



Order No.	Fr	Length (cm)	Guidewire	Qty. / Box
TS6/24	6	24	0.038"	5
TS6/26	6	26	0.038"	5
TS6/28	6	28	0.038"	5
TS7/24	7	24	0.038"	5
TS7/26	7	26	0.038"	5
TS7/28	7	28	0.038"	5

*All products may not be available in all countries. Please contact our sales department for details.

*** Combinations :**

- With Guidewire : Add Suffix - G
- Each stent is supplied with pusher