



allSoft Ureteral Stents (Open Tip, Closed Tip and Intraoperative)



The allSoft Ureteral Stents are intended to temporarily facilitate the drainage of urine from kidney to bladder by a trained physician.

- Low encrustation tendency to facilitate easy removal of the stent
- Soft polyurethane material ensures high patient comfort
- Wide lumen for superior drainage
- Radiopaque stent to enable accurate positioning
- Centimetre marking to aid in stent advancement and precise placement



Order No. (Open Taper Tip)	Order No. (Closed Taper Tip)	Order No. Intraoperative (Both Ends Closed Taper)	Fr	Length (cm)	Guidewire	Qty. / Box
SOT4.7/16	SS4.7/16	SSI4.7/16	4.7	16	0.035"	10
SOT4.7/24	SS4.7/24	SSI4.7/24	4.7	24	0.035"	10
SOT4.7/26	SS4.7/26	SSI4.7/26	4.7	26	0.035"	10
SOT4.7/28	SS4.7/28	SSI4.7/28	4.7	28	0.035"	10
SOT4.7/30	SS4.7/30	SSI4.7/30	4.7	30	0.035"	10
SOT5/16	SS5/16	SSI5/16	5	16	0.035"	10
SOT5/24	SS5/24	SSI5/24	5	24	0.035"	10
SOT5/26	SS5/26	SSI5/26	5	26	0.035"	10
SOT5/28	SS5/28	SSI5/28	5	28	0.035"	10
SOT5/30	SS5/30	SSI5/30	5	30	0.035"	10
SOT6/16	SS6/16	SSI6/16	6	16	0.038"	10
SOT6/24	SS6/24	SSI6/24	6	24	0.038"	10
SOT6/26	SS6/26	SSI6/26	6	26	0.038"	10
SOT6/28	SS6/28	SSI6/28	6	28	0.038"	10
SOT6/30	SS6/30	SSI6/30	6	30	0.038"	10
SOT7/16	SS7/16	SSI7/16	7	16	0.038"	10
SOT7/24	SS7/24	SSI7/24	7	24	0.038"	10
SOT7/26	SS7/26	SSI7/26	7	26	0.038"	10
SOT7/28	SS7/28	SSI7/28	7	28	0.038"	10
SOT7/30	SS7/30	SSI7/30	7	30	0.038"	10

*All products may not be available in all countries. Please contact our sales department for details. Contact us for other sizes

* Combinations :

- With Guidewire : Add Suffix - G
- With Suture : Add Suffix - S
- With Clamp : Add Suffix - C
- With Hydrophilic Coating : Add Suffix - H