


BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip



The BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe
- Bulb shaped tip of the outer sheath
- Silicone depth positioner for precise placement
- The bulb-shaped angled outer sheath eases passage through the cervix



Stylet for BULB TRANS Embryo Transfer Catheter-Bulb Tip

Order No.	Description	Qty. / Box
BT-ETC-U-P	BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip	10
STY-BU	Stylet for BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip	10

*All products may not be available in all countries. Please contact our sales department for details.



Premium bulb tip catheter
with soft inner catheter