



ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush



The ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush is used for laparoscopic or ultrasound - guided transvaginal retrieval of oocytes from ovarian follicles.

- Aspiration and flushing line is made possible through a single channel
- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Dismantling of luer lock connector
- Ergonomic handle
- 100 cm aspiration line




Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
A-OVM-SL16/35-EF	ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush	16	35	10
A-OVM-SL17/35-EF	ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush	17	35	10
A-OVM-SL18/35-EF	ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush	18	35	10

*All products may not be available in all countries. Please contact our sales department for details.



Dismantling of
luer lock connector



Echogenic tip