



ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line




The ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Aspiration and flushing line is made possible through a single channel
- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line
- Multiple operation functionality by using stopcock valve



Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
A-OVM-SL16/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	16	35	10
A-OVM-SL17/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	17	35	10
A-OVM-SL18/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	18	35	10

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip