


ACE-M Manual Single Lumen Ovum Pickup Needle



The ACE-M Manual Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from follicles.

- When used with ultrasonic equipment, echo tip enhances visualization
- Handle designed for greater control during ovum retrieval
- Luer lock connector for connecting a syringe
- The amount of suction applied is manually controlled


Tri-facet geometry
with a echogenic tip



Ergonomic
handle

Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
A-OVM-SL16/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	16	35	10
A-OVM-SL17/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	17	35	10
A-OVM-SL18/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	18	35	10

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip