


## ACE Single Lumen Ovum Pickup Needle



The ACE Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line

Tri-facet geometry  
with a echogenic tip

Ergonomic handle

Aspiration line


Silicone bung  
for tight fit

Luer lock connector

Vacuum line

Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
A-OVM-SL16/35	ACE Single Lumen Ovum Pickup Needle	16	35	10
A-OVM-SL17/35	ACE Single Lumen Ovum Pickup Needle	17	35	10
A-OVM-SL18/35	ACE Single Lumen Ovum Pickup Needle	18	35	10

\*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip