

AspiPump 90D

- High vacuum, low flow pump capable of operating at normal vacuum levels ranging from -10mmHg to -300mmHg and boost vacuum levels from -300mmHg to -500mmHg, with adjustable vacuum settings in increments of -1mmHg.
- Red LED display for easy readability.
- Equipped with a dual foot pedal featuring color coding: Yellow for boost operation and blue for normal operation.
- The Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors (200cm), and Aspi Hydrophobic Filter Set are provided sterile & intended for single use only.

Aspi Suction Bottle with Connector



**COMING
SOON**

ThermiSure TBT-12
(Warmer)

Aspi Hydrophobic Filter Set



Aspi Suction Line with
Luer Lock Connectors



AspiPump Foot Switch



Order No.	Description	Qty. / Box
ASPI 90D	AspiPump 90D Low flow / high vacuum DC operated pump, dual foot pedal, Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors & Aspi Hydrophobic Filter Set, Universal adaptor, User guide Manual.	1
ASP-SB	Aspi Suction Bottle with Connector Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units.	4
ASP-HFS	Aspi Hydrophobic Filter Set Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units.	4
ASP-SL	Aspi Suction Line with Luer Lock Connectors Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units.	4
ASP-FS	AspiPump Foot Switch Suitable for Allwin AspiPump 90D.	4

*All products may not be available in all countries. Please contact our sales department for details.